




# Wearable Sensors in Health Applications

Kamil Grondys


Samsung Electronics, Poland



1. Application Types
2. Web vs Native Sensor API
3. Service Application
4. NPAPI & NPRuntime
5. Gear S2


# Agenda

---



1. Application Types
2. **Web vs Native Sensor API**
3. Service Application
4. NPAPI & NPRuntime
5. Gear S2

# Agenda



1. Application Types
2. Web vs Native Sensor API
- 3. Service Application**
4. NPAPI & NPRuntime
5. Gear S2



1. Application Types
2. Web vs Native Sensor API
3. Service Application
4. **NPAPI & NPRuntime**
5. Gear S2


# Agenda

---




1. Application Types
2. Web vs Native Sensor API
3. Service Application
4. NPAPI & NPRuntime
5. **Gear S2**

# Application Types



# Application Types







1. Magnetic
2. Ambient Light
3. Barometer
4. Proximity
5. Ultraviolet
6. Accelerometer
7. Gyroscope
8. Pedometer
9. Heart Rate Monitor

# Sensor API

---



1. Magnetic
2. Ambient Light
3. Barometer
4. Proximity
5. Ultraviolet
6. Accelerometer
7. Gyroscope
8. Pedometer
9. Heart Rate Monitor

# Sensor API

---



1. Magnetic
2. Ambient Light
3. Barometer
4. Proximity
5. Ultraviolet
6. Accelerometer
7. Gyroscope
8. Pedometer
9. Heart Rate Monitor

# Sensor API

---



1. Magnetic
2. Ambient Light
3. Barometer
4. Proximity
5. Ultraviolet
6. Accelerometer
7. Gyroscope
8. Pedometer
9. Heart Rate Monitor

## How to check which sensors are available on the device:

```
sdb shell -d
```

```
cat /etc/config/model-config.xml
```

```
...  
<key name="tizen.org/feature/sensor.accelerometer" type="bool">true</key>  
<key name="tizen.org/feature/sensor.barometer" type="bool">false</key>  
<key name="tizen.org/feature/sensor.gyroscope" type="bool">true</key>  
<key name="tizen.org/feature/sensor.magnetometer" type="bool">false</key>  
<key name="tizen.org/feature/sensor.photometer" type="bool">false</key>  
<key name="tizen.org/feature/sensor.proximity" type="bool">false</key>  
<key name="tizen.org/feature/sensor.tiltmeter" type="bool">false</key>  
<key name="developer.samsung.com/tizen/feature/heart_rate_monitor" type="bool">true</key>  
<key name="developer.samsung.com/tizen/feature/pedometer" type="bool">true</key>
```


# Sensor API: Example apps

---


1. Fitness measurness (running, walking, etc.)
2. Navigation (distance)
3. Calories intake/outtake

# Sensor API: example for epilepsy management

- Data from the accelerometer



- Epilepsy episode detection



# Sensor API: Tizen API

---

```
enum SensorType {  
    "LIGHT",  
    "MAGNETIC",  
    "PRESSURE",  
    "PROXIMITY",  
    "ULTRAVIOLET",  
    "HRM_RAW"  
};
```

1. MAGNETIC – Magnetic
2. LIGHT – Ambient Light
3. PRESSURE – Barometer
4. PROXIMITY – Proximity
5. ULTRAVIOLET – Ultraviolet
- ...
9. HRM\_RAW – Heart Rate



# Sensor API: HRM\_RAW example

```
var HRMSensor hrm = tizen.sensorservice.getDefaultSensor("HRM_RAW");  
  
hrm.start(function(){  
  
    //on success  
    hrm.getHRMRawSensorData(function(hrmSensorData){  
        //callback to be invoke when sensor data has been read  
        //do something with a data from the sensor  
        //...  
    });  
  
});
```



# Sensor API: HRM\_RAW example

```
var HRMSensor hrm = tizen.sensorservice.getDefaultSensor("HRM_RAW");  
  
hrm.start(function(){  
  
    //on success  
    hrm.getHRMRawSensorData(function(hrmSensorData){  
        //callback to be invoke when sensor data has been read  
        //do something with a data from the sensor  
        //...  
    });  
  
});
```



# Sensor API: HRM\_RAW example

```
var HRMSensor hrm = tizen.sensorservice.getDefaultSensor("HRM_RAW");  
  
hrm.start(function(){  
  
    //on success  
    hrm.getHRMRawSensorData(function(hrmSensorData){  
        //callback to be invoke when sensor data has been read  
        //do something with a data from the sensor  
        //...  
    });  
  
});
```



# Sensor API: Code example for Gear 2



```
var hrmSensor = webapis.motion.getDefaultSensor("HRM");
hrmSensor.start("HRM", function(){
  //on success
  hrmSensor.getHRMRawSensorData(function(hrmSensorData) {
    //sensor data success callback to be invoked when sensor data has been read
    //do something with a data from the sensor
    //...
  });
});
```

## 6. Accelerometer

```
function handleMotionEvent(event)
{
  var x = event.accelerationIncludingGravity.x;
  // or event.acceleration to receive data excluding gravity
  var y = event.accelerationIncludingGravity.y;
  var z = event.accelerationIncludingGravity.z;
  // Do something awesome.
}

window.addEventListener("devicemotion", handleMotionEvent, true);
```

## 7. Gyroscope

```
function handleDeviceOrientationEvent(event)
{
  var alpha = event.alpha;
  var beta = event.beta;
  var gamma = event.gamma;
  // Do something awesome.
}
```

```
window.addEventListener("deviceorientation", handleDeviceOrientationEvent, true);
```

# Sensor API: Pedometer Web vs Native

- Web

```
enum PedometerStepStatus {  
  
    "NOT_MOVING",  
    "WALKING",  
    "RUNNING"  
  
};
```

- Native

```
typedef enum  
{  
  
    ACTIVITY_STATIONARY,    /**< Stationary */  
    ACTIVITY_WALK,         /**< Walking */  
    ACTIVITY_RUN,          /**< Running */  
    ACTIVITY_IN_VEHICLE,   /**< In a moving vehicle */  
  
} activity_type_e;
```

# Sensor API: Native advantages

---

- Possible to get data all the time

```
typedef enum
```

```
{
```

```
    SENSOR_OPTION_DEFAULT,      /**< Does not receive data when the LCD is off  
    */
```

```
    SENSOR_OPTION_IN_SCREEN_OFF, /**< Receives data when the LCD is off */
```

```
    SENSOR_OPTION_IN_POWERSAVE_MODE, /**< Receives data in the PS mode */
```

```
    SENSOR_OPTION_ALWAYS_ON,    /**< Always receives data */
```

```
} sensor_option_e;
```



# Sensor API: Native advantages (set interval)

- Possible to get data all the time

```
typedef enum
```

```
{
```

```
    SENSOR_OPTION_DEFAULT,      /**< Does not receive data when the LCD is off  
    */
```

```
    SENSOR_OPTION_IN_SCREEN_OFF, /**< Receives data when the LCD is off */
```

```
    SENSOR_OPTION_IN_POWERSAVE_MODE, /**< Receives data in the PS mode */
```

```
    SENSOR_OPTION_ALWAYS_ON,    /**< Always receives data */
```

```
} sensor_option_e;
```

- Possible to get even more accurate data

```
int sensor_listener_set_interval(sensor_listener_h listener,  
    unsigned int interval_ms);
```

**Privilege:** <http://tizen.org/privilege/healthinfo>


**Privilege:** <http://tizen.org/privilege/medicalinfo>

# Web vs Native Sensor API: Summary

---

1. JavaScript for Web Sensor API
2. `tizen.sensorservice` since Tizen 2.3
3. `webapis.motion` for older devices
  
4. c/c++ for Native Sensor API
5. Additional options (interval, when to collect data)
6. Higher frequency rate
  
7. Don't forget to add privs

# Web Service Application



1. Runs in the background
2. node.js is used as runtime
3. Partner Level cert. is required
4. Gear 2 not supported

# Web Service Application: when to use

---

1. Data processing without blocking UI
2. Can be launched during a device boot
3. Can communicate with UI application

<https://developer.tizen.org/development/tutorials/web-application/tizen-features/service-application>

# Web Service Application: features

---

```
<widget>
  <tizen:service id="[App_ID]" auto-restart="true" on-boot="true">
    <tizen:content src="[Start_JS_File]"/>
    <tizen:name>[App_Name]</tizen:name>
    <tizen:icon src="[App_Icon]"/>
    <tizen:description>[Description]</tizen:description>
  </tizen:service>
</widget>
```

# Web Service Application: features

---

```
<widget>
  <tizen:service id="[App_ID]" auto-restart="true" on-boot="true">
    <tizen:content src="[Start_JS_File]"/>
    <tizen:name>[App_Name]</tizen:name>
    <tizen:icon src="[App_Icon]"/>
    <tizen:description>[Description]</tizen:description>
  </tizen:service>
</widget>
```

# Web Service Application: example

---

```
module.exports.onStart = function() {
    console.log("service start");
    ...
}

module.exports.onRequest = function() {
    console.log("service received appControl request");
    ...
}

module.exports.onExit = function() {
    console.log("service terminate");
    ...
}
```



# Web Service Application: Summary

---

1. JavaScript service application
2. Can be launched automatically
3. For Partners only
4. Supported since Gear S



1. Plugins for web
2. Available on web since Netscape Navigator web browser
3. c/c++ available for JavaScript
4. Partner Level cert. required

# NPPlugin: how the plugins works

---

- HTML5 code:

```
<object id="plugin" type="application/x-plugin-example"></object>
```

- JavaScript code:

```
var plugin = document.getElementById('plugin')
```

```
plugin.someMethod();  
console.log(plugin.someProperty);
```

# NPPlugin: why?

---

1. Data from a sensor can be obtained with higher freq. rate
2. c/c++ can be used for data processing
3. AI and gesture recognition algorithms can be used
4. Custom codecs can be implemented for Web Apps

# NPPlugin & NPRuntime: Summary

---

1. Call c++ code from JavaScript
2. Tizen Native API available
3. Partner level cert. required
4. Available in Gear 2, Gear S



1. Tizen **Wearable Native API** available
2. Tizen Wearable Web API available
3. Tizen 2.3.1

# Summary

---

1. Health applications often collect data from sensors
2. Sometimes that apps need to work in the background
3. Sometimes an application need to work all the time
4. Wearables can collect and process data
5. You can use both HTML5/JavaScript and C/C++

# Resources

---

1. [developer.tizen.org](http://developer.tizen.org)
2. [developer.samsung.com](http://developer.samsung.com)

 [KamilGrondys](#)

 [KamilGrondys](#)

[k.grondys@samsung.com](mailto:k.grondys@samsung.com)